



WEYERGANS

future health

MICROPROCESSOR-CONTROLLED LIPOLYSIS AND MORE




BODYSTAR

ELECTRO CONTOUR SYSTEM


The whole world of Electro Therapy

Areas of use: for cellulite and local adipose; body contouring and relaxing electro massage of the whole body. For health centres, spas and fitness, cosmetic and wellness institutes.

Technical data: 230 V / 50-60 Hz, capacity 18 W, 8 bit microprocessor, 8 channels, adjustable as a whole and individually via practical remote control, patients' remote control, five different treatment programmes: lipolysis (medical) for needle electrodes, lipolysis (cosmetical) for surface electrodes, polyshape: multi-signal passive gymnastics, lipolysis: local hypertoning and contouring, Alley-Gates: relaxation massage. Can be updated. Dimensions: 430 x 120 x 350 mm (WHD). Robust plastic casing. Colour: blue / grey. Softline-Design as desired. Complete set with accessories and manual.




Polyshape:
A better figure with multi-signal passive gymnastics
Polyshape:
Mit multi-Signal passiver Gymnastik zur besseren Figur



Isotone:
For fast local hypertoning and body contouring
Isotone:
Lokale Hypertonisierung und Body Contouring in kürzester Zeit

As desired Softline side plates in different individual colours available.

Auf Wunsch Softline Seitenteile in verschiedenen individuellen Farben erhältlich.



The wide range of possible applications and uses offers great benefits to the user. Future-oriented state-of-the-art technology due to the possibility of updating the microprocessor! No supervision necessary. Fast amortisation. Lipolysis is a component of differentiated cellulite treatment (DCT).

Bodystar:
Only one device - the whole world of electro therapy

Production and sales:

Weyergans High Care AG
Gut Boisdorf 8 · D-52355 Dueren · Germany
Tel. +49-(0) 24 21-96 78-0
Fax +49-(0) 24 21-96 78-20
info@weyergans.de · www.weyergans.com